



VITROS[®] System
Integrated 5600



A new commitment to innovation

The new OCD's commitment to you, and the patients you serve



- Providing reliable systems with proven technologies that increase productivity and deliver uncompromising result quality.
- Unlocking the future of our proven MicroSlide Technology.
- Re-focus on core menu and investing in novel markers that give clinicians the right knowlege.
- Increasing flexibility and being easy to do business with.
- Seamless integration into your existing lab infustructure.
- Eliminate liquid waste and protect the water supply with systems that don't require water or drains.
- Provide industry-leading uptime, supported by:
 - 1) Proactive monitoring
 - 2) Call centers
 - 3) Experienced field engineers



For nearly 70 years, Ortho Clinical Diagnostics has provided the global healthcare community with the means to make better-informed patient decisions. The OCD of today is re-committed to help your lab deliver the highest quality results to those that need it most.

5 Enabling technologies.


High quality, easy to use, with unmatched flexibility in one system.



VITROS® MicroSlide Technology
Get industry leading 95% reportable result efficiency, with storage and stability ahead of anyone else.



VITROS® MicroWell Technology
Enhanced Chemiluminescence, with exceptional Immunoassay accuracy.



VITROS® MicroTip Technology
Special chemistry menu offerings with less waste, more flexibility.

Know the unknowns.




VITROS® MicroSensor Technology
Minimize interferences and save time, cost, and frustration.

Second guess less.



VITROS® Intellicheck® Technology
Increases result certainty and minimize re-testing.


A better integrated solution has dry slide technology.



**VITROS® is the only solution
with dry slide technology**

No water, plumbing or drains.
You can place it where you want, not where you have to.
Safety in the event of disaster that compromises water supply.

WATER USE BY SYSTEM (🏠 WATER TOWERS/YEAR)



	VITROS® Systems	NON-VITROS® Systems
Water	0 gal/year	43,104 gal/year
Filtration System	None	Required
Plumbing	None	Required
Drains	None	Required
Waste Water Impact to the Environment	None	Present
Additional Cost beyond System	\$0	\$\$\$\$

*Water towers hold 3,500 gallons of water

*CAP Today July, 2013 www.captodayonline.com/productguides

Does your integrated solution
play nice with other systems?



Begin treatment with the right result.



Ortho Clinical Diagnostics Integrated VITROS 5600[®] System delivers the ability to begin treating the world's most serious diseases.

***Our menu covers 90% of the menu needs and 99% of the testing volume for typical labs.**

* Ortho Clinical Diagnostics internal data on file.

VITROS® Menu of Assays

VITROS® ECI/ECiQ Immunodiagnostic Systems				
VITROS® 3600 Immunodiagnostic System				
VITROS® 350 Chemistry System				
VITROS® 4600/5,1 FS Chemistry Systems				
VITROS® 5600 Integrated System				
Cardiac				
Apolipoprotein A1 ²	+	+		
Apolipoprotein B ³	+	+		
Cholesterol ¹	+	+	+	
CKMB	+	+	+	
CKMB Mass	+			+
Direct HDL ¹	+	+	+	
Direct LDL ¹	+	+		
Homocysteine ^{4, 5}	+	+		
hsCRP	+	+		
LDH	+	+	+	
Myoglobin	+			+
Total CK	+	+	+	
NT-proBNP	+			+
Triglycerides	+	+	+	
Troponin I ES	+			+
Drugs of Abuse (Urine)				
Amphetamines	+	+		
Barbiturates	+	+		
Benzodiazepines	+	+		
Cocaine Metabolites	+	+		
DRI® Ecstasy*	+	+		
Methadone	+	+		
Opiate	+	+		
DRI® Oxycodone*	+	+		
PCP	+	+		
THC	+	+		
Infectious Disease				
Anti-HIV 1+2	+			+
Anti-HCV	+			+
Anti-HBc	+		+	+
Anti-HBc IgM	+		+	+
HBsAg	+		+	+
Anti-HBs	+		+	+
Anti-HAV IgM	+		+	+
Anti-HAV Total	+		+	+
Rubella IgG	+		+	+
HBeAg	D		D	+
Anti-HBe	D		D	+
HIV Combo	D		D	D
Syphilis	D		D	D

+ = Available
D = In development/Currently not available
*Manufacturer Validated Assay (MVA); Reagents and validated applications available directly from Thermo Fisher Scientific 1-800-232-3342; product availability subject to third party's or their agent's capability to support the customer.
**Manufacturer Validated Assay (MVA); Reagents and validated applications available directly from Ortho Clinical Diagnostics; product availability subject to third party's or their agent's capability to support the customer.

VITROS® ECI/ECiQ Immunodiagnostic Systems				
VITROS® 3600 Immunodiagnostic System				
VITROS® 350 Chemistry System				
VITROS® 4600/5,1 FS Chemistry Systems				
VITROS® 5600 Integrated System				
Oncology				
AFP	+			+
CA 125 II TM	+			+
CA 15-3 TM	+		+	+
CA 19-9 TM	+		+	+
CEA	+		+	+
IgA ⁷	+	+		
IgG ⁸	+	+		
IgM ⁹	+	+		
PSA	+			+
Total PSA II	D		D	D
Toxicology				
Acetaminophen	+	+	+	
Alcohol	+	+	+	
Cholinesterase	+	+	+	
Salicylate	+	+	+	
Prenatal				
Anti-HIV 1+2	+			+
HBsAg	+			+
Rubella IgG	+			+
Total β-hCG Generation II	+			+
Therapeutic Drug Monitoring				
QMS® Amikacin*	+	+		
Caffeine	+	+		
Carbamazepine	+	+	+	
Digoxin	+	+	+	
Gentamicin	+	+		
Lithium	+	+	+	
EMIT™ Phenobarbital**	+	+		
Phenytoin	+	+	+	
Theophylline	+	+	+	
Tobramycin	+	+		
Valproic Acid	+	+		
Vancomycin	+	+		
Respiratory				
Caffeine	+	+		
Carbon Dioxide	+	+	+	
Chloride	+	+	+	
Potassium	+	+	+	
Sodium	+	+	+	

VITROS® ECI/ECiQ Immunodiagnostic Systems				
VITROS® 3600 Immunodiagnostic System				
VITROS® 350 Chemistry System				
VITROS® 4600/5,1 FS Chemistry Systems				
VITROS® 5600 Integrated System				
Skeletal/Bone				
ALKP	+	+	+	
Calcium	+	+	+	
Intact PTH	+			+
NTx	+			+
Phosphorus	+	+	+	
Vitamin D	+			+
Pancreatic				
Amylase	+	+	+	
Lipase	+	+	+	
Nutritional Assessment				
Albumin ¹⁰	+	+	+	
Direct TIBC	+	+		
Ferritin	+			+
Folate ¹³	+			+
Iron	+	+	+	
Prealbumin	+	+		
Total Protein	+	+	+	
Transferrin ¹²	+	+		
Vitamin B12	+			+
Anemia				
Direct TIBC	+	+		
Ferritin	+			+
Folate ¹³				+
Haptoglobin	D	D		
Iron	+	+	+	
TIBC			+	
Transferrin ¹²	+	+		
Vitamin B12	+			+
General Chemistry				
Amylase	+	+	+	
Calcium	+	+	+	
Glucose	+	+	+	
LDH	+	+	+	
Lipase	+	+	+	
Magnesium	+	+	+	
Phosphorus	+	+	+	
Total CK	+	+	+	
Urea	+	+	+	
Uric Acid	+	+	+	

VITROS® ECI/ECiQ Immunodiagnostic Systems				
VITROS® 3600 Immunodiagnostic System				
VITROS® 350 Chemistry System				
VITROS® 4600/5,1 FS Chemistry Systems				
VITROS® 5600 Integrated System				
Renal				
Albumin ¹⁰	+	+	+	
Calcium	+	+	+	
Carbon Dioxide	+	+	+	
Chloride	+	+	+	
Creatinine	+	+	+	
Glucose	+	+	+	
Microalbumin	+	+		
NephroCheck®	D			D
Phosphorus	+	+	+	
Potassium	+	+	+	
Sodium	+	+	+	
Urea	+	+	+	
Uric Acid	+	+	+	
Urine Protein	+	+	+	
Hepatic				
AAT	+	+		
Albumin ¹⁰	+	+	+	
ALKP ¹¹	+	+	+	
ALT ¹¹	+	+	+	
Ammonia	+	+	+	
AST ¹¹	+	+	+	
BuBc	+	+	+	
Cholinesterase	+	+	+	
GGT	+	+	+	
IgA	+	+		
IgG	+	+		
IgM	+	+		
Total Bilirubin	+	+	+	
Prealbumin	+	+		
Spinal				
Glucose (csf)	+	+	+	
Total Protein (csf)	+	+	+	
Lipids				
Apolipoprotein A1	+	+		
Apolipoprotein B	+	+		
Cholesterol	+	+	+	
Direct HDL	+	+	+	
Direct LDL	+	+		
Triglycerides	+	+	+	

VITROS® ECI/ECiQ Immunodiagnostic Systems				
VITROS® 3600 Immunodiagnostic System				
VITROS® 350 Chemistry System				
VITROS® 4600/5,1 FS Chemistry Systems				
VITROS® 5600 Integrated System				
Inflammatory				
ASO ¹⁵	+	+		
C3	+	+		
C4	+	+		
CRP ¹⁴	+	+	+	
IgA	+	+		
IgG	+	+		
IgM	+	+		
RF	+	+		
Thyroid/Metabolic				
Cortisol	+			+
Free T3	+			+
Free T4	+			+
Intact PTH	+			+
Lactate	+	+	+	
T3 Uptake	+			+
Total T3	+			+
Total T4	+			+
TSH 3 rd Generation	+			+
Diabetes				
Creatinine	+	+	+	
Glucose	+	+	+	
HbA1c	+	+		
Microalbumin ⁶	+	+		
Urea	+	+	+	
Reproductive				
Estradiol	+			+
FSH	+			+
LH	+			+
Progesterone	+			+
Prolactin	+			+
Testosterone	+			+
Total β-hCG Generation II	+			+
TSH 3 rd Generation	+			+
Immunosuppressant Drugs				
Cedia® Cyclosporine*	+	+		
QMS® Everolimus*	D	D		
Cedia® Mycophenolic Acid*	+	+		
QMS® Tacrolimus*	D	D		

References

1. Warnick GR, Myers GL, Cooper GR, Rifai N. Impact of the Third Cholesterol Report from the Adult Treatment Panel of the National Cholesterol Education Program on the clinical laboratory. *Clin Chem.* 2002;48(1):11-17.

2. Ritchie RF (ed). *Serum Proteins in Clinical Medicine*, Vol 1. Scarborough, ME: Foundation for Blood Research; 1996:12.01-12.05.

3. Sandkamp M. Apolipoprotein-Diagnostik: klinische Relevanz. (in German). *Diagnose & Labor.* 1990;40:37-43.

4. Sundström J, Sullivan L, D'Agostino RB, et al. Plasma homocysteine, hypertension incidence, and blood pressure tracking: the Framingham Study. *Hypertension.* 2003;42(6):1100-1105.

5. Zylberstein DE, Bengtsson C, Bjorkelund C, et al. Serum homocysteine in relation to mortality and morbidity from coronary heart disease: a 24-year follow-up of the population study of women in Gothenburg. *Circulation.* 2004;109(5):601-606.

6. Burtis C, Ashwood ER (eds). *Tietz Textbook of Clinical Chemistry.* 3rd ed. Philadelphia, PA: WB Saunders Co; 1995:484.

7. Tietz NW, Pruden E, McPherson RA, Fuhrman SA (eds). *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia, PA: WB Saunders Co; 1995:355-357.

8. Tietz NW, Pruden E, McPherson RA, Fuhrman SA (eds). *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia, PA: WB Saunders Co; 1995:358-361.

9. National Committee on Clinical Laboratory Standards. Protection of laboratory workers from instrument biohazards and infectious diseases transmitted by blood, body fluids, and tissue: approved guideline. Wayne, PA: National Committee on Clinical Laboratory Standards; 1987.

10. Tietz NW (ed). *Fundamentals of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders Co; 1987:328-329.

11. Tietz NW (ed). *Fundamentals of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders Co; 1987:369-371.

12. Burtis C, Ashwood ER (eds). *Tietz Textbook of Clinical Chemistry.* 3rd ed. Philadelphia, PA: WB Saunders Co; 1995:500-501.

13. Chanarin I. *The Megaloblastic Anemias.* 3rd ed. Oxford, United Kingdom: Blackwell Science; 1990:6.

14. Gambino R. C-reactive protein—undervalued, underutilized. *Clin Chem.* 1997;43(11):2017-2018.

15. Ayoub EM, Harden E. Immune response to streptococcal antigens: diagnostic methods. In: *Manual of Clinical Laboratory Immunology*, 6th ed. Rose NR, Hamilton RG, Detrick B, eds. Washington, DC: American Society of Microbiology; 2002.

NephroCheck® and the NephroCheck® logo are registered trademarks of Astute Medical, Inc. in the United States. CA 125 II, CA 15-3, and CA 19-9 are trademarks of Fujirebio Diagnostics, Inc. CEDIA®, DRI® and QMS® are registered trademarks of Thermo Fisher Scientific and its subsidiaries. HIV and Hepatitis co-developed with Chiron/Novartis. www.orthoclinical.com All trademarks are the property of their respective owners. © Ortho-Clinical Diagnostics, Inc. 2009-2015 PR-00224

Ortho Clinical Diagnostics

Integration that is truly integrated.



Saves time and money

Chemistry and immunoassay samples are processed in parallel. Improves turn around time and increases efficiency.



More checks less balances

Matched processing capabilities for chemistry and immunoassay with uninterrupted workflow.



More checks less samples

150 reagent positions and 100 assays onboard at one time.



Don't Stop.

Matched processing capabilities for chemistry and immunoassay with uninterrupted workflow.



No carry-over EVER

Single use tips for sample and reagent metering.



Don't call us we'll call you

e-Connectivity® / Proactive Monitoring

